

根据不同的样品量，NORELL 制作出适用于 Bruker MATCH™ 系统的八种超精密核磁管，根据管帽的颜色可以辨别出管子的型号，为用户提供便利。

产品编号	外径 (mm)	内径 (mm)	长度 (mm)	样品体积 (μl)	包装数量 (支)
S-1.0-500-1	1.00+0.010-0.025	0.58 ± 0.010	100	12	5
S-1-0.73-500-1	1.00+0.010-0.025	0.73 ± 0.010	100	19	5
S-1.7-500-1	1.70+0.010-0.025	1.30 ± 0.010	100	45	5
S-2.0-500-1	2.00+0.010-0.025	1.60 ± 0.010	100	70	5
S-2.5-500-1	2.50+0.010-0.025	2.10 ± 0.010	100	120	5
S-3.0-500-1	2.99+0.010-0.025	2.41 ± 0.010	100	160	5
S-4.0-500-1	3.99+0.010-0.025	3.20 ± 0.010	100	310	5
S-4.25-500-1	4.25+0.010-0.025	3.43 ± 0.010	100	370	5
S-5.0-500-1	4.97+0.010-0.025	4.20 ± 0.010	100	490	5

